BUNSEN

USER'S MANUAL

BATHS FOR OIL SAMPLING

BATHS FOR OIL SAMPLING

This equipment is designed for sensory analysis of olive oil by the method "Organoleptic assessment of virgin olive oil" (COI/T.20/Doc. n°15/Rev.1)

This method has been reviewed and approved by the International Olive Council in accordance with standard COI/T.20/Doc.n°13/Rev.1

The bath is formed by:

- Digital thermostatic bath
- Aluminum block 45 x 240 x 240 mm

TECHNICAL FEATURES

- Working temperature; adjustable by means of digital thermostat from room temperature +5° to 99.9°C.
- Temperature reading by means of a 3 digit display.
- Resolution; 0,1°C
- Accuracy; ±0,5% to end of scale
- Temperature control; ON/OFF mode
- Probe type PTC
- Capacity 9 litres
- Power 1200 W

CHARACTERISTICS

- Power supply voltage; 230v-50Hz.
- Equipment designed for indoor use only.
- Double body metal construction, anti-rust treated exterior, epoxy painted and stainless steel 18/10 vat.
- Heating systems by means of AISI 316 stainless steel plated elements.
- Equipped with security thermostat covering excess temperature cases on the rear side of the equipment (DIN 12877 CLASS 2).
- Power supply cord, Black PVC (approximate diameter 7.2 mm) H05VV-F3G - 1.0 mm cordage meets CENELEC HD-21 (HAR) requirements. 2 meters long.
- European "SCHUKO" angle plug with 4,8 mm pins. 90 angle cord grip. Meets CEE 7-7 requirements, DIN 49440 / 41. "SCHUKO" outlets and connectors (EU1-16R).

INSTALLATION

The equipment should be placed on a flat, levelled and sturdy surface.

The supply connection is by means of a 2m long power supply cable with side earth connection plug with 4.8mmØ "SCHUKO" type locking bananas.

IMPORTANT!

In order to avoid the possible appearance of rust on the bottom of the vat, please clean this with paper or a slightly moist cloth before filling the equipment with water.

SAFETY RECOMMENDATIONS

Follow the safety recommendations to prevent damage to persons or property.

- Connect the unit only to a grounded mains power socket.
- Operation is permitted with non-flammable liquids only.
- Never operate the unit without bath fluid in the bath.
- Be careful when emptying hot bath fluids. Check the temperature of the bath fluid before draining.
- Never operate damaged or leaking equipment.
- Always turn off the unit and disconnect the mains cable from the power source before performing any service or maintenance procedures or before moving the unit.
- Always empty the bath before moving the unit.
- Never operate equipment with damaged mains power cables.
- Some parts of the cover bath may become extremely warm during continuous operation. Therefore, exercise particular caution when touching these parts.

CLEANING THE UNIT


Before cleaning the unit, disconnect the power plug from the mains socket.

Sometimes it's happens that if water of the network is used, depositions of lime in the bath tank and in the resistance be produced. You must avoid that the resistance work cover with lime, because the resistance heats less and its useful life diminishes.

To clean of lime the bath tank to fill it of water with a 10% rate of nitric acid, until all the lime be removed. After empty it and rinse it with clean water.

FUNCTIONING

Control panel

1. Heating indicating light (amber coloured).
 2. Safety thermostat operation indicating light.
 3. Start light switch.
 4. Digital thermostat.
- For your safety, do not start the bath if this is not connected to a good earth connection.
- Before connecting the equipment to the electricity supply, make sure that the voltage on the characteristics plate is the same as that of the electric power supply.
- Fill the bath with liquid, making sure than its level is always above the equipment element.
- Connect the equipment to the power supply.
- Put the start-up switch into position **I**.

- Programme the digital thermostat at the temperature you desired to work. (See digital thermostat user manual).
- During the liquid heating-up period of time, the amber-coloured light will be switched on indicating that the element is working. Once the working temperature is reached, this will be switched off, switching on and off intermittently in order to maintain said temperature stable.
- Take into account the liquid evaporation; try not to leave the bath element uncovered.

PROGRAMMING OF THE THERMOSTAT

1. DESCRIPTION OF THE CONTROL PANEL

- 1. Set button:** Allows the selection of the Set point and the operation parameters.
- 2. DOWN ▼ button:** Allows the selection of said programming parameters and decrease their value.
- 3. UP ▲ key:** Allows the selection of the various programming parameters and the increase in their value.
- 4. F key:** Programmable key not used in this version.
- 5. Led Set:** Indicates that the programming modality and the parameter programming level are activated.
- 6. Led Out:** Indicates the thermostat exit status (functioning of the element); switched on (heats up); switched off (does not heat up).
- 7. Led Def:** Not enabled.

The thermostat is configured by default in such a way that the user has only to choose the temperature to work on.

2. PROGRAMMING OF THE SET POINT

Press button “**Set**” and the letters `sp` will be shown in the display (in the factory this value is 10, in order to prevent the element from starting to heat up when accidentally starting the equipment) alternating with the set value.

In order to modify the programmed value press key \blacktriangle or \blacktriangledown depending on whether the value to be worked at is desired to be increased or decreased.

These keys act digit by digit, although if they remain pressed for more than a second, the value increases or decreases quickly, and after being pressed for two seconds, the speed increases in order to quickly reach the desired value.

Once the desired value is programmed, by pressing the “**Set**” key you exit the fast programming mode.

In order to exit the Set fast programming mode, press key “**Set**” after the visualisation of the last Set, or, either do not press any key for 15 seconds, after which the display will return to its normal operation mode.

3. PARAMETER PROGRAMMING

In this section how to access all the configuration parameters of the thermostat is explained, although there are a series of parameters that should not be modified by the user, as this can provoke a malfunctioning of the equipment, although there are some, such as the `ofs` calibration parameter, which may be interesting so the user can check the temperature measured by this equipment with a reference thermometer the user may have. Access to these parameters is carried out in the following way.

- Keep the **Set** button pressed for at least 5 seconds until the value `0` is shown.
- Press button `▲` until value `381` is shown, which is the configuration access password.
- Press button **Set** and parameter `spl1` will be shown.
In order to change the parameter value, press button **Set** and a numeric value alternating with `spl1` will be shown (in this case). Then press `▲` or `▼` to modify the value of the parameter.
- Once we have chosen the desired value, by pressing **Set** again only the parameter `SPL1` will be shown.
- In order to continue visualising the following parameters press button `▲` or `▼` until visualising the one to be modified and follow the same procedure as that when modifying the `spl1` value.

The order of appearance of the parameters, starting with spl1 and pressing button ▼, is the following:

1. spl1
2. sph1
3. sens
4. ofs
5. dp
6. hset
7. func
8. od
9. PAss
10. SP

There are a series of parameters, specially designed for cooling, which is not shown as the thermostat is only programmed to work as a heater.

TABLE OF PARAMETERS

PARAMETERS	DESCRIPTION	RANGE	VALUE
1 SPLL	Minimum Set point	-504SPHL	0
2 SPHL	Maximum Set point	SPLL4300	99,9
3 SenS	Type of probe	Ptc-ntc	Ptc
4 OFS	Probe calibration	-30430 °C	0
5 dP	Decimal point	On-Off	On
6 HSEt	Differential	0,0430,0°C	0,0
7 Func	OUT functioning mode HEAT-COLD	Heat-Cool	Heat
8 PASS	Access password to operation parameters	OFF49999	381
9 SP	Set Point	SPLL4SPHL	10,0

ERROR SIGNALS

ERROR	CAUSE	SOLUTION
E1-E1	The probe is stopped or short circuited or a value outside the allowed range has been measured	Check the correct connection of the probe with the instrument and the correct functioning of the probe
EEPr	Internal memory error	Verify, and if necessary, re-programme the operation parameters

4. PARAMETERS PROGRAMMING USING THE PASSWORD

The instrument has a parameter protection function using a password that can be personalised through the "pass" parameter

If one wishes to have this protection (pass = xxxx, 4 numbers from 0000 to 9999) set the password number desired in the parameter "pass" following the next steps:

- Press the **SET** key to access the parameters and keep it press for about 5 seconds, after which the LED SET will flash and the display will show 0.
- Push the ▲ or ▼ keys until set the password number
- Push the **Set** key and the display will show spl1.
- To change the parameter value push the **Set** key and will show a number alternating with spl1.
- Push the ▲ or ▼ keys to modify the value.
- When the value is chosen push the **Set** key again and will show the parameter spl1.
- If you want to visualise other parameters push the ▲ or ▼ keys until show the value that you want modify. Then the parameter can be modified using the same programming procedure to program spl1.
- The parameters will show in the same way described in the previous section.
- Protection using a password can be disabled by setting the parameter "pass"="off".

SAFETY THERMOSTAT

The device includes a safety thermostat, adjustable from 100° to 320°C, on the rear side of the bath.

Its use is to protect the element in the event of overheating, normally due to a lack of liquid in the bath. In the event this is activated, the red light will be switched on and the element will be disconnected avoiding the operation of the bath (does not heat up). In order to start the bath again, you have to wait for the temperature to decrease by around 50%; then press the re-arm switch on said thermostat (small red button).

Depending on the maximum working temperature of the bath (100° or 200°C) the safety thermostat will be thus regulated by default, the user being able to change the same, if so desired, by following the scale marked on the thermostat, always taking into account that an excessive safety temperature may seriously damage the bath element.

GUARANTEE

This equipment has a **2 YEAR** guarantee covering spare parts and labour costs for any manufacturing defect, as from the date of purchase written down by the seller.

The breakdowns caused by inadequate use, handling errors or any other reason not attributable to a manufacturing defect are excluded from this guarantee.

For any query or request for spare parts, the equipment model and the series number has to be provided (located on the characteristics plate).

BUNSEN, S. A. CANNOT BE HELD RESPONSIBLE FOR ANY DAMAGE THAT MAY RESULT FROM INADEQUATE USE, HANDLING OR OPERATION OF THE EQUIPMENT COVERED BY THIS GUARANTEE.

Product manufactured by:

BUNSEN, S. A.

Tel. 34916113584

Fax. 34916128254

S.A.T: mam@bunsen.es

EC DECLARARTION OF CONFORMITY

MANUFACTURER: BUNSEN, S. A.

ADDRESS: Tokio, 2. Pol. Ind Wells. 28970 Humanes de Madrid (Madrid – Spain)

PRODUCT:

WATER BATH BD SERIES. Ref. 1625

**TO WHICH THIS DECLARATION RELATES IS IN CONFORMITY WITH THE
FOLLOWING STANDARS:**

EN 61010-2-010

EN 61000-4-3

EN 55014.

EN 61000-3-2

EN 61000-3-3

STATES:

That the product described above fulfil the European Directives:

Directive for machinery 73/23/EEC with amendment 93/68/EEC.

Directive electromagnetic compatibility 89/336/EEC with amendments 91/263/EEC,
93/31/EEC and 93/68/EEC.

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