cryokit

CRYOGENIC PROTECTION

NOTA INFORMATIVA E ISTRUZIONI PER L'USO
NOTE D'INFORMATION ET INSTRUCTIONS D'UTILISATIONS
NOTA INFORMATIVA E INSTRUCCIONES DE USO
INFORMATION NOTE AND INSTRUCTIONS FOR USE
INFORMATIONSBLATT UND GEBRAUCHSANWEISUNG

CRYOPLUS-2.0
CRYOPLUS-2.1
CRYOKIT 300-400-550
CRYOLITE-WP-HP-HPS

RELEASE 1/2022

INFO@CRYOKIT.NET

CRYOKIT.NET

CRYOPLUS-2.0
CRYOPLUS-2.1

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CRYOPLUS-2.0
CRYOPLUS-2.1

CRYOKIT300
CRYOKIT400
CRYOKIT550

These gloves are made with a special composite fabric bonded to a hydrophilic membrane. The interior features another full, breathable membrane, which does not allow liquids to penetrate inside. Protection from contact cold is guaranteed by an internal lining in soft polyester fleece.

CE CERTIFICATION

CE certification issued by: Ricotest Srl – Notified Lab No. 0498
Subject to specific liquid nitrogen resistance test (-195.82°Cc), by immersion in contact for 1 minute and subsequent RCT test.

N.B.: the performance level concerns the palm side, including fingers; nevertheless, the uniformity of materials and processing makes the back of the glove also protective.

The following marking can be found on the label sewn into the inside of the gloves:

These gloves are made with special composite fabrics bonded to hydrophilic membranes. The interior features another special, breathable membrane, which does not allow liquids to penetrate inside. Protection from contact cold is guaranteed by an internal lining in soft polyester fleece.

CE CERTIFICATION

CE certification issued by: Ricotest Srl – Notified Lab No. 0498
Subject to specific liquid nitrogen resistance test (-195.82°Cc), by immersion in contact for 1 minute and subsequent RCT test.

N.B.: the performance level concerns the palm side, including fingers; nevertheless, the uniformity of materials and processing makes the back of the glove also protective.

The following marking can be found on the label sewn into the inside of the gloves:

	CE marking, testifying to the conformity of the gloves to the essential health and safety requirements of Regulation (EU) 2016/425
0302	Identification no. of the Notified Body performing the inspection of the PPE manufactured pursuant (MODEL C2)
CRYOKIT	Manufacturer's registered trade mark
CRYOPLUS-2.0/2.1	Glove model
10 (example)	Size
	Pictogram referring to the information note.
abcde	EN388+pictogram relating to the mechanical risks explained below: 2(a) 3(b) 4(c) 2(d) x(e) This pictogram indicates that the gloves were designed for mechanical risks and the numbers at the foot indicate the levels obtained during laboratory tests (see explanation on the last page).
abc	EN 511 + related pictogram for protection against cold explained below: (a) 2 = level obtained against convective cold (0.193 m ² K/W) (b) 2 = level obtained against contact cold (0.1563 m ² K/W) (c) 1 = level obtained against water penetration (no penetration)

	CE marking, testifying to the conformity of the gloves to the essential health and safety requirements of Regulation (EU) 2016/425
0498	Identification no. of the Notified Body performing the inspection of the PPE manufactured pursuant (MODEL C2)
CRYOKIT	Manufacturer's registered trade mark
CRYOKIT400-550	Glove model
10 (example)	Size
	Pictogram referring to the information note.
abcde	EN388+pictogram relating to the mechanical risks explained below: 2(a) 2(b) 3(c) 1(d) x(e) This pictogram indicates that the gloves were designed for mechanical risks and the numbers at the foot indicate the levels obtained during laboratory tests (see explanation on the last page).
abc	EN 511 + related pictogram for protection against cold explained below: (a) 2 = level obtained against convective cold (0.157 m ² K/W) (b) 2 = level obtained against contact cold (0.0562 m ² K/W) (c) 1 = level obtained against water penetration (no penetration)

CRYOLITE-WP
CRYOLITE-HP

These gloves are made of soft water-repellent cowhide grain leather. The interior features another special, breathable membrane, which does not allow liquids to penetrate inside. Protection from contact cold is guaranteed by an internal lining in soft polyester fleece. **CE CERTIFICATION**
CE certification issued by: Ricotest Srl – Notified Lab No. 0498
Subject to specific liquid nitrogen resistance test (-195.82°Cc), by immersion in contact for 1 minute and subsequent RCT test.

N.B.: the performance level concerns the palm side, including fingers; nevertheless, the uniformity of materials and processing makes the back of the glove also protective.
The following marking can be found on the label sewn into the inside of the gloves:

	CE marking, testifying to the conformity of the gloves to the essential health and safety requirements of Regulation (EU) 2016/425
0498	Identification no. of the Notified Body performing the inspection of the PPE manufactured pursuant (MODEL C2)
CRYOKIT	Manufacturer's registered trade mark
CRYOLITE	Glove model
10 (example)	Size
	Pictogram referring to the information note.
 abcde	EN388+pictogram relating to the mechanical risks explained below: 2(a) 1(b) 2(c) 2(d) x(e) This pictogram indicates that the gloves were designed for mechanical risks and the numbers at the foot indicate the levels obtained during laboratory tests (see explanation on the last page).
 abc	EN 511 + related pictogram for protection against cold explained below: (a) 1 = level obtained against convective cold (0.224 m ² K/W) (b) 3 = level obtained against contact cold (0.1476 m ² K/W) (c) 1 = level obtained against water penetration (no penetration)
	EN 407 + related pictogram for protection against heat explained below: X2XXXX 2 = Level obtained against contact heat: 16 seconds at 250°C.

CRYOLITE
HPS

These gloves are made of soft, yet tough, water-repellent cowhide grain leather. The interior features another special, breathable membrane, which does not allow liquids to penetrate inside. Protection from contact cold is guaranteed by an internal lining in soft polyester fleece. **CE CERTIFICATION**
CE certification No. 617162701/OI issued by: Ricotest Srl – Notified Lab No. 0498
Subject to specific liquid nitrogen resistance test (-195.82°Cc), by immersion in contact for 1 minute and subsequent RCT test.

N.B.: the performance level concerns the palm side, including fingers; nevertheless, the uniformity of materials and processing makes the back of the glove also protective.
The following marking can be found on the label sewn into the inside of the gloves:

	CE marking, testifying to the conformity of the gloves to the essential health and safety requirements of Regulation (EU) 2016/425
0498	Identification no. of the Notified Body performing the inspection of the PPE manufactured pursuant (MODEL C2)
CRYOKIT	Manufacturer's registered trade mark
CRYOLITE HPS	Glove model
10 (example)	Size
	Pictogram referring to the information note.
 abcde	EN388+pictogram relating to the mechanical risks explained below: 3(a)32(b) 3(c) 2(d) x(e) This pictogram indicates that the gloves were designed for mechanical risks and the numbers at the foot indicate the levels obtained during laboratory tests (see explanation on the last page).
 abc	EN 511 + related pictogram for protection against cold explained below: (a) 2 = level obtained against convective cold (0.216 m ² K/W) (b) 2 = level obtained against contact cold (0.1517 m ² K/W) (c) 1 = level obtained against water penetration (no penetration)
	EN 407 + related pictogram for protection against heat explained below: X2XXXX 2 = Level obtained against contact heat: 16 seconds at 250°C.

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Kora Srl – Cryokit Cryogenic