Peltier vs Compresseur

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The background of the slide features a close-up photograph of a leaf's vein structure. A semi-transparent green rectangular overlay covers the lower half of the image, with its right edge rounded. The word 'eco' is printed in white on this green area. The Memmert logo is positioned in the top right corner.

memmert

trust the best

eco


Efficient. Homogeneous. Reliable

HPP_{eco} / IPP_{eco} | PELTIER TECHNOLOGY

Memmert has been the expert in the application of Peltier technology in laboratory equipment for over 20 years.

Maximum energy efficiency
compared to climate chambers
with a compressor.

Absolutely homogeneous
distribution of temperature
and humidity.

Reliable and fail-safe,
even after many years
of continuous operation.

Always a good climate:
refrigerant-free, smooth running,
compact, low maintenance.

HPP_{eco} / IPP_{eco} | ADVANCED PELTIER TECHNOLOGY

New Peltier technology

- Highly efficient in continuous operation
 - Twice the power with the same energy consumption
 - High savings on Total Cost of Ownership
-


New micro-metering pump

- Intelligent humidity control
 - Fast, accurate and reproducible
-

New air flow routeing

- Outstanding homogeneity
due to modified air flow routeing
-

HPPeco / IPPeco | ENERGY SAVING

HPP_{eco} / IPP_{eco} | EXAMPLE CALCULATION - ICH GUIDELINE

Energy consumption for long-term tests (12 months) according to ICH-Q1A (R2)

Energy costs

Energy costs

HPP410_{eco}

Compressor device

Calculation:

438 KWh x € 0.25

5200 KWh x € 0.25

Energy costs:

€ 109.5

€ 1300

Environmental Impact

CO₂ calculation

HPP410_{eco}

Compressor device

Calculation:

438KWh x 0,537kgCO₂ (KWh)

5200KWh x 0,537kgCO₂ (KWh)

CO₂ Footprint:

0,24 t CO₂

2,8 t CO₂

HPP^{eco} / IPP^{eco} | ADVANCED PELTIER TECHNOLOGY

New Peltier technology

- Highly efficient in continuous operation
 - Twice the power with the same energy consumption
 - High savings on Total Cost of Ownership
-

New micro-metering pump

- Intelligent humidity control
 - Fast, accurate and reproducible
-

New air flow routeing

- Outstanding homogeneity
due to modified air flow routeing
-

HPP_{eco} / IPP_{eco} | NEW HUMIDIFICATION SYSTEM

New vaporisation plant

- Enables faster humidification



New micro-metering pump

- Enables 10 times more precise humidification

HPP110_{eco}*

Spatial humidity deviation	25 °C/ 60% rh	+ - 0.2 %rh
Temporal humidity deviation	25 °C/ 60% rh	+ - 0.7 %rh
Humidification at 35 °C	10% rh to 60% rh	8 min

* The values are based on measurements on a pre production series. Subject to modifications.

HPP^{eco} / IPP^{eco} | ADVANCED PELTIER TECHNOLOGY

New Peltier technology

- Highly efficient in continuous operation
 - Twice the power with the same energy consumption
 - High savings on Total Cost of Ownership
-

New micro-metering pump

- Intelligent humidity control
 - Fast, accurate and reproducible
-

New air flow routeing

- Outstanding homogeneity
due to modified air flow routeing
-

HPPeco / IPPeco | NEW AIR FLOW ROUTEING


- Outstanding homogeneity due to changed air flow routeing in the interior
- Faster recovery times after opening the door
- Optimal temperature distribution, even at maximum load

HPPeco / IPPeco | APPLICATIONS

- Long term stability tests
- Packaging tests
- Quality checks
- Humidity-controlled storage

” **Ideal for stability tests
according to ICH Q1A**

HPP^{eco} / IPP^{eco} | TECHNICAL DATA

Model sizes: 110, 260, 410, 750, 1060, 1400, 2200

Temperature adjustment ranges: 0 to +70 °C (110 -1060) | +15°C to +60°C (1400/2200)

Humidity adjustment ranges: 10 to 90% rh (HPP110^{eco} - HPP1060^{eco}) | 10 to 80% rh (HPP1400^{eco}/HPP2200^{eco})