

memmert

trust the best



Water bath

WTB

WTB | PRACTICAL AND ROBUST

Easy to clean

- rounded corners
- Stainless steel surfaces

Easy to empty

- new drain system

Easy to handle

- smallest possible footprint




WTB | INTUITIVE AND USER-FRIENDLY

Easy to operate

- 3.5 inch touchscreen
- New user interface


Versatile functions

- Timer
- Alarm
- 2-point calibration
- Setting Celsius / Fahrenheit



WTB | CUSTOMVIEW

Thanks to the customisable remote View **CustomView** you always keep an eye on the important information.



WTB | TECHNICAL DATA



WTB6
7.5 l



WTB11
10 l



WTB15
15 l



WTB24
24 l



WTB35
37.5 l



WTB50
51 l

Model sizes / Article description	6	11	15	24	35	50
Power consumption at 230/115 V, 50/60 Hz (W)	1000		1420		2800	
Operating temperature range (°C)	min. 5 above ambient temperature up to +100					
Setting temperature range (°C)	+10 to +100					